

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922) 49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Ижевск (3412)26-03-58
Иваново (4932)77-34-06
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Сургут (3462)77-98-35
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольяти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Улан-Удэ (3012)59-97-51
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://cleaver.nt-rt.ru> || cea@nt-rt.ru

Центрифуги



Quickspin Mini Centrifuge

SKU: QUICKSPINPLUS

- Supplied with both standard microtube and strip tube rotors
- Ideal for quick spin-downs and microfiltration

- Starts and stops in seconds
- Compact design

The Quickspin is perfect for microfiltration and rapid spin-down of sample from the walls and caps of microcentrifuge tubes.

By occupying less than 6 inches square of bench space, the Quickspin has a very small footprint, making it easy to use in the lab. Rotors and adaptors, which may be easily interchanged, are supplied as standard to accommodate 1.5 ml, 0.5 ml and 0.4 ml tubes, as well as 0.2 ml strips and individual tubes.

A highly durable stainless steel hinge pin facilitates easy opening of the translucent lid, while an on/off switch is located on the side of the centrifuge to start and stop operation. Alternatively, with the switch in the on position, the centrifuge can be started and stopped by closing and opening the lid.

SPECIFICATIONS

Unit Dimensions (cm)	15 x 15 x 11.6
Weight (Kg)	1.2
Speed	6000 rpm
Capacity	8 x 1.5/2.0ml, 4 x 0.2ml strips
Includes Rotor	Yes
Maximum G	2000
Deceleration Time	<2 seconds
Input Voltage	100-240V



multiFUGE Mini Centrifuge

SKU: CSL-MULTIFUGE

Unlike traditional mini centrifuges, the new MultiFUGE Mini Centrifuge eliminates the need to change rotors when switching between microtubes and PCR strips; the unique DuoROTOR is all that is required for running 12 microtubes and 4 PCR strips simultaneously.

With a fixed speed that produces 2,000 xg, the MultiFUGE is perfect for quick spin downs. Simply close the lid and the unit quickly ramps up to 5500 rpm. Open the lid and the rotor quickly decelerates for removal of samples.

At just 14 cm wide and 11 cm high, the MultiFUGE Mini Centrifuge is a compact centrifuge with unmatched capacity and flexibility.



QUICKSPIN2 Mini-Centrifuge

SKU: **QUICKSPIN2**

- A locking lid for added protection when the unit is in use
- A larger lid to accommodate longer and skirted tubes, as well as taller tube lids
- A braking function to quickly stop the unit when powered off or when the lid is opened
- Meets the latest applicable North American and EU safety requirements

SPECIFICATIONS

Maximum speed: Approx. 6,000 rpm (fully loaded)

Maximum RCF: 2000 X g

Power: 100 to 240 VAC 50/60 HZ 0.6A, Approx. 7.3 W

multiFUGE Adaptors, 0.2ml, Pack of 6

SKU: **MF-A0.2-6**

multiFUGE Adaptors, 0.5ml, Pack of 6

SKU: **MF-A0.5-6**



- **ezeSPIN Adaptor, 0.2ml PCR tube (8/pk)**

SKU: EZ-A0.2



- **ezeSPIN Adaptor, 0.5ml PCR tube (8/pk)**

SKU: EZ-A0.5



•

ezeSPIN Microcentrifuge Tube Rotor, 8 x 1.5/2.0ml

SKU: EZ-ROTOR



•

ezeSPIN mini centrifuge with digital speed adjustment

SKU: EZEESPIN-D

- Quick release rotor system – no tools needed to exchange rotors
- Electronic braking – Allows for more rapid deceleration, less waiting
- Higher speed rating than comparable models
- Maximum capacity: 8 x 1.5/2.0 ml (standard rotor); 4 x strip tubes (strip tube rotor)
- Rubber coated base plate offers improved stability and reduced vibration
- ezeSPIN Digital comes with digital display and control for both speed and timer

SPECIFICATIONS

EZEEspin-D

Speed (rpm)	1,000 - 12,000 rpm
Max RCF	6,209 × g
Max Capacity	8 × 2ml
Set Time	1 to 99mins
Power Supply	90-240VAC, 50/60Hz
Includes	2 Rotors, 0.5ml & 0.2ml Adaptors
Weight, Kg	1.328
Dimensions (WxDxH) mm	176x196x122



- **ezeSPIN PCR Strip Tube Rotor, 4 x 8 PCR strip**

SKU: **EZ-PST**



- **ezeSPIN single speed mini centrifuge**

SKU: **EZEESPIN**

- [SPECIFICATIONS](#)
-

Speed (rpm)	7,000 rpm
Max RCF	2,113 × g
Max Capacity	8 × 2ml
Set Time	NA
Power Supply	90-240VAC, 50/60Hz
Includes	2 Rotors, 0.5ml & 0.2ml Adaptors
Weight, Kg	1.324
Dimension (WxDxH) mm	176x196x122



•
High Speed Centrifuge and a Magnetic Stirrer In-One

SKU: MAGFUGE

- Magnetic Stirrer and High Speed Centrifuge All-In-One
- Low noise level helps users stay focused in the lab
- Long-life, maintenance free brushless motor
- High quality rare-earth magnets in stirrer rotor ensure strong magnetic coupling
- Two magnetic stir bars are included to help ensure maximum stirring efficiency



24 place microcentrifuge (Includes 24×1.5/2.0ml rotor)

SKU: EZEEFUGE

- Large touchscreen display provides clear indication of running conditions.
- Up to 20,625 x g, tube capacity up to 12 x 5 ml.
- Supports labware including PCR tubes, strip tubes, 2ml microtubes, 5ml tubes.
- Store up to 40 programs with a combination of speed, centrifugal force, running time, acceleration & deceleration rate.
- 3 available rotors



24 place refrigerated microcentrifuge

SKU: See Product Page

- Large touchscreen display provides clear indication of running conditions.
- Up to 20,625 x g, tube capacity up to 12 x 5 ml.
- Supports labware including PCR tubes, strip tubes, 2ml microtubes, 5ml tubes.
- Store up to 40 programs with a combination of speed, centrifugal force, running time, acceleration & deceleration rate.
- 3 available rotors
- Temp range -20°C~40°C



Fixed Angle Rotor 24 x 1.5/2ml, 15,000 rpm

Алматы (7273)495-231
 Ангарск (3955)60-70-56
 Архангельск (8182)63-90-72
 Астрахань (8512)99-46-04
 Барнаул (3852)73-04-60
 Белгород (4722)40-23-64
 Благовещенск (4162)22-76-07
 Брянск (4832)59-03-52
 Владивосток (423)249-28-31
 Владикавказ (8672)28-90-48
 Владимир (4922) 49-43-18
 Волгоград (844)278-03-48
 Вологда (8172)26-41-59
 Воронеж (473)204-51-73
 Екатеринбург (343)384-55-89

Ижевск (3412)26-03-58
 Иваново (4932)77-34-06
 Иркутск (395)279-98-46
 Казань (843)206-01-48
 Калининград (4012)72-03-81
 Калуга (4842)92-23-67
 Кемерово (3842)65-04-62
 Киров (8332)68-02-04
 Коломна (4966)23-41-49
 Кострома (4942)77-07-48
 Краснодар (861)203-40-90
 Красноярск (391)204-63-61
 Курск (4712)77-13-04
 Курган (3522)50-90-47
 Липецк (4742)52-20-81

Киргизия (996)312-96-26-47

Магнитогорск (3519)55-03-13
 Москва (495)268-04-70
 Мурманск (8152)59-64-93
 Набережные Челны (8552)20-53-41
 Нижний Новгород (831)429-06-12
 Новокузнецк (3843)20-46-81
 Ноябрьск (3496)41-32-12
 Новосибирск (383)227-86-73
 Омск (3812)21-46-40
 Орел (4862)44-53-42
 Оренбург (3532)37-68-04
 Пенза (8412)22-31-16
 Петрозаводск (8142)55-98-37
 Псков (8112)59-10-37
 Пермь (342)205-81-47

Россия (495)268-04-70

Ростов-на-Дону (863)308-18-15
 Рязань (4912)46-61-64
 Самара (846)206-03-16
 Саранск (8342)22-96-24
 Санкт-Петербург (812)309-46-40
 Саратов (845)249-38-78
 Севастополь (8692)22-31-93
 Симферополь (3652)67-13-56
 Смоленск (4812)29-41-54
 Сочи (862)225-72-31
 Ставрополь (8652)20-65-13
 Сыктывкар (8212)25-95-17
 Сургут (3462)77-98-35
 Тамбов (4752)50-40-97
 Тверь (4822)63-31-35

Казахстан (772)734-952-31

Тольяти (8482)63-91-07
 Томск (3822)98-41-53
 Тула (4872)33-79-87
 Тюмень (3452)66-21-18
 Улан-Удэ (3012)59-97-51
 Ульяновск (8422)24-23-59
 Уфа (347)229-48-12
 Хабаровск (4212)92-98-04
 Чебоксары (8352)28-53-07
 Челябинск (351)202-03-61
 Череповец (8202)49-02-64
 Чита (3022)38-34-83
 Якутск (4112)23-90-97
 Ярославль (4852)69-52-93